



## Minncare® Cold Sterilant **Research Data**



**MAR COR**  
**PURIFICATION**  
A Cantel Medical Company



# Minnicare® Cold Sterilant Research Data

## CONTENTS

### EFFECTIVENESS

In Vitro Aqueous Testing .....4

AOAC Sporicidal Testing .....5

Vapor Kill Testing .....6

### TOXICITY ASSESSMENT SUMMARY

Acute Oral Toxicity Testing .....7

Inhalation Toxicity Testing .....7

Testing of Skin Sensitivity .....8

Mucous Membrane Sensitivity Tests .....9

Dermatological Sensitivity Tests .....9

Environmental Effects .....10

### STABILITY

Stability Over Time .....11

Temperature Stability .....12

Storage Conditions .....13

MATERIAL SAFETY DATA SHEETS .....14

## EFFECTIVENESS

### In Vitro Aqueous Testing

To show that Minncare is effective against other microorganisms, the following in vitro tests were performed using a .5% solution of Minncare. All tests were done on aqueous test solutions at 20°C.

#### In Vitro Aqueous Test at 20°C

Minncare Concentration .5%		
Species	Count per ml	Time for 100% kill in minutes
Bac. subtilis	$6 \times 10^6$	2.5
Bac: stearothermophilus	$6 \times 10^6$	2.5
Bac. subtilis NCTC 3610	$2.4 \times 10^9$	5.0
Bac. mesentericus	$1.6 \times 10^9$	5.0
Clostr. perfringens	$1 \times 10^7$	10.0
Clastr. tyrobutyricum	$1 \times 10^7$	5.0
Sacchar. cereisiae	$6 \times 10^7$	0.5
Cand. mycoderma	$1.4 \times 10^8$	0.5
Hansenula anomala	$6.4 \times 10^8$	0.5
Pichia membronaefaciens	$4.8 \times 10^8$	0.5
Pen. camerunense	$1.7 \times 10^8$	2.5
Mucor plumbeus	$3 \times 10^6$	2.5
Geotrichum candidum	$2 \times 10^7$	0.5
Byssoschlamys nivea	$6 \times 10^7$	0.5
Staph. aureus	$2.6 \times 10^9$	0.5
Strept. faecalis	$4.6 \times 10^9$	0.5
Kleb. aerogenes	$2.3 \times 10^9$	0.5
Ps. fluorescens	$4.6 \times 10^9$	0.5
Ps. aeruginosa	$2 \times 10^9$	0.5
Salm. thyphimurium	$2.8 \times 10^9$	0.5
Coryneb. rubrum	$1 \times 10^7$	1.0
Leuconostoc spec.	$5.3 \times 10^8$	0.5
Lactob. brevis	$1.8 \times 10^9$	0.5

## AOAC Sporidical Testing

To show that Minncare is sporidical, the test method of the American Organization of Analytical Chemist (AOAC) was performed. This test is accepted by the E.P.A. to demonstrate the efficacy of a chemical as a sterilant.

Cultures of *Bacillus subtilis* and *Clostridium sporogenes* were grown. Both silk suture and loops and ceramic cylinders were contaminated with the cultures and dried for 24 hours in a vacuum oven. Tests were then run to show that after dehydration the spores are still viable and that resistance is within tolerance.

Five loops or cylinders were then placed into 10ml of the test solution, in this case diluted Minncare solution. The test time was set for 11 hours.

Following the 11 hour contact time, the carriers were removed from the test solution and individually placed into test tubes containing a subculture medium. The test samples were then transferred to a fresh tube of thioglycolate and incubated for 21 days at 37°C.

A total of 720 tests were performed on three lots of the test chemical. Results of the test are shown on the chart below.

AOAC Sporidical Testing		
Lot/Organism	Carrier Type	Positives/Total
1 <i>Bacillus subtilis</i>	loop	0/60
1 <i>Bacillus subtilis</i>	cylinder	0/60
1 <i>Clostridium sporogenes</i>	loop	0/60
1 <i>Clostridium sporogenes</i>	cylinder	0/60
2 <i>Bacillus subtilis</i>	loop	0/60
2 <i>Bacillus subtilis</i>	cylinder	0/60
2 <i>Clostridium sporogenes</i>	loop	0/60
2 <i>Clostridium sporogenes</i>	cylinder	0/60
3 <i>Bacillus subtilis</i>	loop	0/60
3 <i>Bacillus subtilis</i>	cylinder	0/60
3 <i>Clostridium sporogenes</i>	loop	0/60
3 <i>Clostridium sporogenes</i>	cylinder	0/60

### NO POSITIVE CULTURES WERE OBTAINED, THUS MINNCARE PASSED THE AOAC SPORIDICAL TEST.

The performance standards are as follows:

- For sporidical claims, no more than 2 failures can be tolerated.
- For sterilizing claims, no failures can be tolerated.
- Growth must be observed in tubes with carriers exposed for 2 minutes to 2.5N HC (per the dehydration/ resistance test).
- Controls must show growth.

All of these performance standards were met by an 11 hour exposure to Minncare solution.

## Vapor Kill Testing

The ability of Minncare vapor to kill *Bacillus subtilis* spores was tested in vitro.

An aqueous solution of *Bisubtilis* spores was made by grinding spore strips with sterile water in a sterile container. The solution was then filtered through .45 µm membrane filters. A specific amount of the solution was filtered through each membrane to give a 10<sup>3</sup>, 10<sup>4</sup>, 10<sup>5</sup>, 10<sup>6</sup> spores respectively. Each of the filters supported 1" above a 1% solution of Minncare. Controls had the filters supported 1" above sterile water. All jars were then covered lightly with plastic caps. The jars were allowed to stand for 21 hours at ambient temperature.

The filters were then removed and transferred to petri dishes containing plate count agar. The dishes were incubated for 48 hours at 35°C. The results are recorded in the chart below.

Deactivation of Organisms By Minncare Vapor	
Estimated spore concentration of filter suspended 1 inch above	B. subtilis spores/filter
Control: water	loop
10 <sup>3</sup> spores	9.6 x 10 <sup>3</sup> CFU*
10 <sup>6</sup> spores	> 3 x 10 <sup>3</sup> CFU
Minncare solution (1%)	
10 <sup>3</sup> spores	<1 CFU
10 <sup>4</sup> spores	<1 CFU
10 <sup>5</sup> spores	<1 CFU
10 <sup>6</sup> spores	<1 CFU
* Colony forming unit	

## TOXICITY ASSESSMENT SUMMARY

AOAC Sporidical Testing			
Toxicity Summary	Specimen	Results	Reference
Oral toxicity	Male Rats	LD <sub>50</sub> = 2.43 (2.04-2.88) g/kg	Litchfield & Wilcoxon
	Female Rats	LD <sub>50</sub> = 2.10 (1.92-2.30) g/kg	
Inhalation toxicity	Male Rats	Established lethal concentration	
	Female Rats	LC = 13,439 mg/cubic meter	
Skin sensitivity	Mice	No reaction	Burkhard's Test
	Humans	No visiable effects	
Mucos membranes	New Zealand Rabbits	Effects completely gone within 7 days	HH Draiz
Dermatological sensitivity qualities	White Guinea Pigs	No difference between control and test group	Klugman & Magnusson
Intravenous toxicity	Male Rats	LD <sub>50</sub> = 212 mg/kg	Litchfield & Wilcoxon

### Acute Oral Toxicity Testing

The test for acute oral toxicity of Minncare solution was conducted utilizing male wistar rats weighing 197 grams and female wistar rats weighing 157 grams. Ten rats were used per dosage. The Minncare was diluted with water and administered to the animals by means of force feeding hose. The solution was administered at a constant 20 ml/kg of body weight. The animals were observed after the treatment for period of 8 days. As symptoms of poisoning, the following variables were observed: skin rashes, decreased mobility, difficulties in breathing, cramps, and inability to stand up. The LD50 value was statistically, it was shown to be:

Male rats: LD<sub>50</sub> = 2.43 (2.04 - 2.88) g/kg

Female rats: LD<sub>50</sub> = 2.10 (1.92 - 2.30) g/kg

The dosage that was survived by all rats was 1.25 g/kg for male rats and 1.58 g/kg for female rats. The difference is due to the weight difference between males and females. After testing was completed, all animals were dissected. An analysis revealed etching of internal organs, but indicated that there were no specific toxic effects.

### Inhalation Toxicity Testing

Because the possibility exists that persons exposed to Minncare for even a short period of time could inhale vapors, we have tested Minncare solution for acute inhalation toxicity with animals. Ten male and ten female rats were placed in an inhalation room and given a 5% solution of Minncare in the form of a fine mist. This concentration was inhaled by the rats for a period of 4 continuous hours without any toxic symptoms.

---

The test was then repeated with another 20 rats. This time, undiluted Minncare was sprayed into the inhalation chamber. Over a 4 hour test period, 38.2 grams of diluted solution and 28.8 grams of undiluted solution were administered to the animals. These products were mixed with 2123 and 2023 liters of air, respectively. The capacity of the inhalation chamber was 120 liters; therefore, 1.91 grams of undiluted Minncare from the diluted solution and 28.8 grams of undiluted Minncare were sprayed into the test chamber during the two tests. This resulted in concentrations of 851 mg/m<sup>3</sup> of Minncare in the test using 5% (diluted) concentration of Minncare and 13,429 mg/m<sup>3</sup> of Minncare in the test using undiluted.

Neither group of rats that were exposed to the Minncare solution showed any symptoms of poisoning. The rats that were sprayed with the undiluted product exhibited the following symptoms: scratching their noses, inducement to sneeze, wet skin, general discomfort (which was expressed by crawling together in a corner of the room), and bent backs. These symptoms remained for 1 hour after the test ended. All 40 rats survived the chamber for the observed period of 1 week.

The lethal concentration of Minncare is therefore established at greater than 13,439 milligrams per cubic meter.

We can conclude from these tests that Minncare in the vapor form, when used according to the instructions, does not create any health hazard for personnel.

## Testing of Skin Sensitivity

### Repeated Application on the Skin of Hairless Mice

A group of 10 hairless mice was tested for skin sensitivity to Minncare. Minncare concentration of 2% was applied to a skin patch the size of a silver dollar twice a day for 2 weeks. It was rubbed into the skin and left there. Two mice showed a slight reddening of the skin after the fourth treatment. After the sixth treatment, this was still the case with 1 mouse. These were the only animals who displayed any such symptoms. The other animals did not show any symptoms during or after the treatment. After the eighth and through the twelfth treatment, none of the mice exhibited any skin reaction.

### Burckhardt - Test with Volunteers (Humans)

The testing of skin sensitivity with repeated applications to 5 human volunteers was conducted according to the W. Burckhardt-Test as described in the professional Dermatology Magazine (p 179-188, 1970). Using a glass stick, a 3% Minncare solution was applied to a skin patch the size of a silver dollar located on the inside of each test subject's underarm. The treatment was repeated at 30 second intervals for a period of 30 minutes.

Neither during nor after this test, were there any visible effects of the product on the skin



## Mucous Membranes Sensitivity Test

Two groups of four white male New Zealand Rabbits were used. Both groups were given 0.1 ml of a 3% water solution of Minncare. The solution was applied to the tear duct of the right eye. The left eye of the rabbits was left untreated to serve as a control. Group one's treated eyes were thoroughly rinsed 10 seconds after the Minncare solution was applied. Group two's eyes were not rinsed after the application.

Per the time schedule developed by the Draize method, the assessment of the effect on the cornea, iris, and the white of the eye followed 2, 6, 24, 48, 72 and 144 hours after the treatment (H.H. Draize, "Appraisal of the Safety of Chemicals in Foods, Drugs and Cosmetics," Assn. of Food and Drug Officials of the U.S., pp. 49-52, 1959). No effects on the cornea or iris and no infectious symptoms were observed. A severe reddening of the white of the eye occurred with both groups. The reddening experienced by groups one and two was 57.5% and 77.5% of the maximum possible reaction, respectively. This reddening decreased greatly within 3 days after the treatment, and was completely gone 7 days after the treatment.

This research indicates that the contact with mucous membranes in the eye with both undiluted or slightly diluted product must be avoided. Contact of the eye with diluted Minncare can result in reddening of the eye. If such contact does occur, a fast, intensive rinsing of the eye should be performed, and an ophthalmologist should be consulted. Safety glasses should always be worn when working with Minncare solution

## Dermatological Sensitivity Test

This test was conducted according to the method of Kligman and Magnusson (Journal of Investigation Dermatology, Volume 52, pp. 268-276, 1979). Two groups of pure white guinea pigs were used. Each guinea pig was between 300-400 grams. A spot 5 x 6 cm was carefully shaved behind their shoulder blades.

The test animals in group one were given three intracutaneous injections on each side of the back bone. The injections were as follows: 1) 0.1 ml Freund' Schem adjuvans placebo; 2) 0.1 ml solution of 1% Minncare; and 3) .1ml of mixture of .05 ml placebo and .05 ml 1% Minncare solution.

The control group two were given two intracutaneous injections on each side of the backbone. The injections were as follows: 1) 0.1 ml placebo; and 2) 0.1 ml distilled water. Eight days later, a 1 % solution of Minncare and Vaseline® was applied to the shaved patch with gauze and held in place with tape. This gauze was removed 48 hours later.

After a 14 day interval in treatment, an area 2 x 2 cm was shaved on the right flank of all two groups of animals. Vaseline® mixed with a 1% Minncare solution was applied to the area and covered with gauze. The gauze was taped in place for 24 hours. After the removal of the gauze, all animals in both group one and group two showed equally strong reddening of the skin, which must be given a Kligman and Magnusson value of 1. Twenty-four hours later, this was only the case with 4 animals. Another 24 hours later, none of the 20 animals showed any reddening of the skin.

---

Because we cannot observe a difference between the test animals and the control animals, it was concluded that Minncare does not have any sensitivity effect on the skin. In this test, only primary skin reactions which were caused by the 1% Minncare solution were observed.

## Environmental Effects


The effects of Minncare on the environment were evaluated according to standard testing methods. (USEPA Guideline No. 72-1) The tests included the following:

1. Minncare as  $H_2O_2$  : Trimetric Analytical Method Validation in Freshwater. Laboratory project #J9207003d.
2. Minncare Acute Toxicity to Rainbow trout, *Oncorhynchus mykiss*, Under Static Test Conditions. Laboratory project #J9207003b.
3. Minncare: Acute Toxicity to Bluegill, *Lepomis macrochirus*, Under Static Test Conditions. Laboratory project #J9207003c.
4. Minncare: Acute Toxicity to the Water Flea, *Daphnia magna*, Under Static Test Conditions. Laboratory project #J9207003a.

The results of this testing are available upon request


## Stability Over Time

Although Minncare will remain stable in its concentrated form for over 1 year, once it is diluted, a decay process begins to take place. The decay is such that when diluted with AAMI quality water, 50% of the peracetic acid will remain after 7 days. The graph that follows shows the decay curve of Minncare diluted 1:100 over a 14 day period. The ratio of dilution effects the rate of decay of peracetic acid contained in the Minncare solution. Thus, if the dilution ratio is less than 1:100, the rate of decay will be less.



## Temperature Stability

Concentrated Minncare remains stable at temperatures up to 30°C (86°F) for a period of 1 year. However, once Minncare is diluted, the rate of decay of the peracetic acid is greatly increased as the temperature increases.



%Peracetic Acid in Diluted Solution

### Minncare Stability:

Once the product is diluted, it must be used within 7 days, and storage and use temperature must be maintained below 30°C (86°F).

## Storage Conditions

A test was conducted on a .5% solution of Minncare to demonstrate the stability of peracetic acid in various storage conditions. The storage conditions researched consisted of translucent containers at room temperature, dark containers at room temperature and dark containers at 37°C. The selection of these test conditions was based on the knowledge that both temperature and light can affect the level of peracetic acid. The results of this test are shown in the table below:

Stability of Peracetic Acid with Storage Conditions

Minncare Concentration .5%			
Storage Condition	% of Decrease in 7 Days	% of Decrease in 14 Days	% of Decrease in 112 Days
Translucent Container Room Temperature	43%	58%	83%
Dark Container Room Temperature	36%	43%	80%
Dark Container 37°C	51%	78%	78%

## Section 1. Product and Company Identification

---

### **Identification:**

Product Name: MINNCARE® COLD STERILANT (EPA Reg. No. 52252-4)

MINNCARE® Liquid Disinfectant (as sold in Canada)

DIN 02277484

Company Identification: Minntech Corporation

14605 28<sup>th</sup> Avenue North

Minneapolis, MN 55447

1-800-328-3340

(763) 553-3300

Emergency Telephone Number: CHEMTREC 1-800-424-9300, or (703) 527-3887

## Section 2. Composition / Information on Ingredients

---

Ingredient	CAS #	Amount (percentage by Weight)	PEL
Hydrogen Peroxide	7722-84-1	22.0%	1 ppm
Peracetic Acid	79-21-0	4.5%	NE
Acetic Acid	64-19-7	-	10 ppm
Water	7732-18-5	-	NE

PEL(s) represent the OSHA 29 CFR 1910.1000 the eight hour time weighted average (TWA) for Hydrogen Peroxide and Acetic Acid.

NE = None Established

## Section 3. Hazards Identification

---

Appearance: Clear

Physical State: Liquid

Odor: Acid

Hazards of Product: Corrosive and an Oxidizer

### **Potential Health Effects**

**Inhalation:** Effect from inhalation of mist will vary from mild irritation to serious damage of the upper respiratory tract, depending on severity of exposure. Symptoms may include sneezing, sore throat or runny nose.

**Ingestion:** Swallowing may cause severe burns of mouth, throat, and stomach. Severe scarring of tissue may result. Symptoms may include bleeding, vomiting, diarrhea, fall in blood pressure.

**Skin Contact:** Contact with skin can cause irritation or severe burns.

**Eye Contact:** Causes irritation of eyes, and with greater exposures it can cause burns that may result in permanent impairment of vision.

## Section 4. First Aid Measures

---

**Eyes and Skin:** Flush with large amounts of water for at least 15-20 minutes. Remove contact lenses, if present, after the first five (5) minutes and then continue rinsing. If burn or irritation has occurred, seek medical attention. If clothing is contaminated, remove clothing, wash skin and wash clothing before reusing.

**Ingestion:** If swallowed, sip a glass of water if able to swallow. Do not attempt to induce vomiting.

**Inhalation:** If inhaled, move to fresh air.

## Section 5. Fire Fighting Measures

---

Flash Point:	N/A
Flammable Limits:	N/A
Extinguishing Media:	Water, Foam CO <sub>2</sub> , Dry Chemicals
Unusual Fire and Explosion Hazards:	N/A

## Section 6. Accidental Release Measures

---

Put on eye protection, protective gloves, boots, clothing and a respirator if air contamination is above the permitted levels. Contain the spill and neutralize with sodium bicarbonate or sodium carbonate. If allowed by federal, state or local regulatory authority, flush spill to the sewer. If mops, towels, paper towel or similar material is used, insure that these items are thoroughly rinsed with copious amounts of water. Do not reuse the liquid material.

## Section 7. Handling and Storage

---

**General Handling:** Keep container closed, but vented when not in use. Store in a cool, dry area (below 75°F). Store unused product in original closed container. Once the product has been removed, do not return to the original container.

## Section 8. Exposure Controls / Personal Protection

---

<b>Eyewear:</b>	ANSI approved safety glasses or goggles. A face shield should be worn when splashes are likely.
<b>Gloves:</b>	Protective gloves should be worn.
<b>Clothing:</b>	A protective apron should be worn when splashes are likely. Rubber boots should be used for spill response.
<b>Respirator:</b>	If air contamination is above the permitted levels, use an NIOSH approved respirator.

## Section 9. Physical and Chemical Properties

---

<b>Physical State:</b>	Liquid
<b>Appearance:</b>	Clear
<b>pH (as a concentrate):</b>	0.5 – 1.1
<b>Solubility in Water (by weight):</b>	Complete
<b>Odor:</b>	Acid
<b>Molecular Weight:</b>	ND
<b>Boiling Point (760 mmHg):</b>	ND
<b>Freezing Point:</b>	ND
<b>Specific Gravity (H<sub>2</sub>O = 1):</b>	1.13
<b>Vapor at Pressure at 20°C:</b>	ND
<b>Vapor Density (air = 1):</b>	ND
<b>Evaporation Rate (Butyl Acetate=1):</b>	ND
<b>Melting Point:</b>	ND

## Section 10. Stability and Reactivity

---

<b>Conditions to Avoid:</b>	Hot storage
<b>Incompatible Materials:</b>	Metals including iron, copper, copper alloys, brass and aluminum, salts, flammable organics, alkalis, caustics, chlorine and formaldehyde.
<b>Hazardous Polymerization:</b>	Will not occur.
<b>Hazardous Decomposition:</b>	Do not mix chlorinated products as this could liberate toxic corrosive chlorine gas.

## Section 11. Toxicological Information

Minnicare® as a product. LC<sub>50</sub> for inhalation is > 2.26 mg/L. LD<sub>50</sub> for oral ingestion is 2.10 g/kg.

### Cancer Related Information:

Ingredient	CAS #	NTP	IARC	OSHA
Hydrogen Peroxide	7722-84-1	Known: No Anticipated: No	None	NO
Peracetic Acid	79-21-0	Known: No Anticipated: No	None	NO
Acetic Acid	64-19-7	Known: No Anticipated: No	None	NO
Water	7732-18-5	Known: No Anticipated: No	None	NO

## Section 12. Ecological Information

**Environmental Fate:** No information found.

**Environmental Toxicity:** This product is toxic to birds, fish and aquatic invertebrates

## Section 13. Disposal Considerations

Dispose of this product in accordance with all applicable, Federal, State and Local regulations.

## Section 14. Transport Information

NON-Bulk

**Proper Shipping Name:** Hydrogen Peroxide and Peroxyacetic Acid Mixtures, Stabilized

**Hazard Class:** Oxidizer (5.1) and Corrosive (8)

**UN Number:** 3149

**Packing Group:** II

## Section 15. Regulatory Information

### International Inventory Status:

Ingredient	CAS #	EC	Japan	Australia	Korea	Canada: DSL	Canada: NDSL
Hydrogen Peroxide	7722-84-1	YES	YES	YES	YES	YES	NO
Peracetic Acid	79-21-0	YES	YES	YES	YES	YES	NO
Acetic Acid	64-19-7	YES	YES	YES	YES	YES	NO
Water	7732-18-5	YES	YES	YES	YES	YES	NO

### United States:

Ingredient	CAS #	OSHA	CAA	CWA	RCRA	SARA 302	SARA 313	TSCA
Hydrogen Peroxide	7722-84-1	YES	NO	NO	NO	NO	NO	NO
Peracetic Acid	79-21-0	YES	YES	NO	NO	YES	YES	NO
Acetic Acid	64-19-7	YES	NO	YES	NO	NO	NO	NO
Water	7732-18-5	YES	NO	NO	NO	NO	NO	NO

**CA Proposition 65:** This product is not affected by CA Proposition 65.

**WHMIS (Canada):** This MSDS has been prepared according to the hazard criteria of the Controlled Products Regulations (CPR) and the MSDS contains all of the information required by CPR.



## Section 16. Other Information

---

### NFPA Ratings:

Flammability: 0  
Health: 2  
Reactivity: 1  
Specific Hazard: Corrosive

### HMIS Ratings:

Flammability: 0  
Health: 2  
Reactivity: 1  
PPE: B

Origination Date: 1/10/90  
Revision Date: 5/13/11  
Prepared by: Corporate Director of Risk Management







Mar Cor Purification  
4450 Township Line Road  
Skippack, PA 19474-1429  
Tel: (484) 991-0220  
Toll Free: (800) 346-0365  
Fax: (484) 991-0230

Mar Cor Purification  
14550 28th Avenue North  
Plymouth, MN 55447  
Tel: (484) 991-0220  
Toll Free: (800) 633-3080  
Fax: (763) 210-3868

Mar Cor Purification  
3250 Harvester Road - Unit 6  
Burlington, ON L7N 3W9  
Tel: (905) 639-7025  
Toll Free: (800) 268-5035  
Fax: (905) 639-0425

Mar Cor Purification  
Sourethweg 11  
6422 PC Heerlen  
The Netherlands  
Tel: (+31) 45 5471 471  
Fax: (+31) 45 5429 695

Mar Cor Purification  
1A International Business  
Park, #05-01  
Singapore 609933  
Tel: (+65) 6227 9698  
Fax: (+65) 6225 6848